



PRODUCT DATA SHEET

RBP 1500A™ Aerosol Activator

Description

- Instantly cures cyanoacrylate adhesive
- Fine atomization and consistent coverage promote more even curing
- The can will not spill or leak like pump type activators
- oz. net wt.

Directions

1. Apply RhinoBond adhesive per instructions
2. Lightly mist RBP 1500A onto the adhesive
3. Be careful not to spray too much activator. This will cause the adhesive to "boil," therefore reducing the bond strength
4. The adhesive should be cured or hardened after 10 seconds
5. reapply with a light mist to cure any uncured adhesive
6. Once adhesive is cured, cut mesh and remove the frame.